

Hi-Neck™

Disk Care


Cervical Traction Device



What is the Hi-Neck™?

The special design of the Hi-Neck™ (Cervical traction device) gives anatomically correct traction. This cervical traction device supports the natural curve of the cervical spine (C-curve).

Also, Hi-Neck™ gently stretches neck muscles allowing the vertebral discs to realign, freeing the nerve root tissue from the pressure of the discs. Hi-Neck™ helps “locked” neck muscles to slowly stretch and relax, allowing misaligned vertebrae to resume their normal supportive position.



Indications

Hi-Neck™ cervical traction pillow helps with conditions below.

- Neck pain
- Whiplash
- Headache
- Occipital neuralgia
- Herniated disc
- Bulging disc
- Degeneration in the neck spine
- Loss of curve in the neck spine
- Chronic neck pain
- Degeneration or arthritis in the neck

What makes this 'Hi-Neck™' a standout?

- User friendly without restraints
- Portable, comfortable and convenient
- Require no assembly, no weights, no cables, no springs or levers. It's simply restraint free
- Encourage regular use by patients at home, helping in breaking the chronic pain cycle. The special design of the Hi-Neck™ gives anatomically correct traction
- Support the natural curve of the cervical spine (C-curve) (Ex. Military neck or Loss of curve in the neck spine)
- Help to solve sleeping disorder without taking medications
- Alleviate neck and shoulder pain resulting from sitting at a computer, being on your feet all day or minor strains in your neck
- Support better circulation because blood will flow better when your muscles are relaxed

Instruction For Use

1. Place the Hi-Neck™ on the floor.
2. Lower into the device ensuring that the bottom of your head can be supported between the contoured lobes.
3. Settle back the full weight of your head, shoulders and neck on the pillow to experience complete comfort and support from the gentle traction force it generates.
4. Apply it for 20~30minutes per time. (a couple of times per day)



Hi-Neck™

Disk Care

Cervical Traction Device