


















FORA-B PLUS

Model information & Spec

No.	Model Name	Catheter		Guide wire	Needle	Port	Syringe
1	SWCB38	SW-BS-03		Length 560 O.D. 1.1			
			SW-BC-32			SWCV-12	HJL3
		COV33B-1		O.D. 2.1			
			SWV33B-1		SWGW-700-085	SWVP-33L	
2	SWCB39	SW-BS-03		Length 560 O.D. 1.1			
			SW-BC-32			SWCV-12	HJL3
		COV33B-2		O.D. 2.1			
			SWV33B-1			SWVP-33L	
3	SWCB40	SW-BS-03		Length 560 O.D. 1.1			
			SW-BC-32			SWCV-12	HJL3
		COV33B-3		O.D. 2.1			
			SWV33B-2		SWGW-700-085	SWVP-33L	




FORA-B PLUS

Minimally invasive
Quick recover
Easy & Safe



FORA-B PLUS can be performed by injecting a drug into the epidural space to relieve pain.
As a one-piece medical device, It is a balloon catheter with a balloon at the end of the catheter.
In addition, the built-in guidewire can be combined to improve the rigidity when inserting the catheter.



FORA-B PLUS



Advantage

- Low risk therapy
- Minimally invasive
- No general anesthetic
- High success rates
- No scar formation
- No open surgery
- No long hospitalization
- Quick recovery
- Repeatable at any time
- Also ideal for patients who have undergone previous surgery
- Targeted treatment of causes

Feature

- Perform bidirectionally
- Radiopaque & soft aromatic Shaft
- Steerable Aromatic distal tip
- Easy to handle
- Easy to steer precisely
- No need to withdrawn the guidewire to inject the medication

Contraindications

- Pregnancy
- Infection
- Mental illness