







RHINO






Disposable Handpieces

This product is a disposable device used to cut soft tissue or bones during surgery. These products are disposable and used in electric medical handpieces.

Burs

Tip / Handle			Length(mm)	O.D
	SWST-BB45-130	SWLV-BB45-130	130mm	Ø4.5
BARREL BURS (5 Flute)	For rapid, aggressive subacromial decompression, notchplasty, and distal clavicle resection without clogging.			
	SWST-RB45-130	SWLV-RB45-130	130mm	Ø4.5
ROUND BURS (8 Flute)	For rapid, aggressive general abrasion arthroplasty, notchplasty, and distal clavicle resection without clogging.			

Cutters

Tip / Handle			Length(mm)	O.D
	SWST-CB42-130	SWLV-CB42-130	130mm	Ø4.2
RESECTOR	For general debridement. Particularly effective on articular cartilage. Ideal for removal of periosteum from the acromion. Run in forward for bone removal in acromioplasty and notchplasty.			
	SWST-JB42-130	SWLV-JB42-130	130mm	Ø4.2
TOMCAT	The narrow gap between inner and outer cutting surfaces provides aggressive resection of soft tissue in knee and shoulder. Extremely effective on meniscus and plicae.			
	SWST-AG42-130	SWLV-AG42-130	130mm	Ø4.2
AGGRESSIVE PLUS	Razor sharp toothed inner blade provides aggressive resection of meniscus, cartilage, and synovium. No teeth on outer window to protect articular cartilage upon entry.			